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1. Shion, H.; *et al.* Enabling Routine and Reproducible Intact Analysis When Data Integrity Matters. Waters Application Note. 2019. 720006472EN <<https://www.waters.com/nextgen/us/en/library/application-notes/2019/enabling-routine-and-reproducible-intact-mass-analysis.html>> .
2. Ranbaduge, N.; *et al.* Routine Peptide Mapping Analysis using the BioAccord System. Waters Technology Brief. 2019. 720006466EN. <<https://www.waters.com/nextgen/us/en/library/application-notes/2019/routine-peptide-mapping-analysis-using-bioaccord.html>>
3. Zhang, X.; *et al.* Released N-linked Glycan Analysis using the BioAccord System. Waters Technology Brief. 2019. 720006474EN. <<https://www.waters.com/nextgen/us/en/library/application-notes/2019/released-n-linked-glycan-analysis-using-bioaccord.html>>
4. Shion, H.; *et al.* A Platform Method for the Molecular Mass Analysis of the Light Chains and Heavy Chains of Monoclonal Antibodies using the BioAccord System. Waters Technology Brief. 2019. 720006529EN. <<https://www.waters.com/nextgen/us/en/library/application-notes/2019/platform-method-molecular-mass-analysis-light-heavy-chains-monoclonal-antibodies-bioaccord-system.html>>
5. Zhang, X.; *et al.* Increasing Productivity and Confidence for N-linked Glycan Analysis of Biosimilars using the BioAccord System. Waters Technology Brief. 2019. 720006545EN. <<https://www.waters.com/nextgen/us/en/library/application-notes/2019/increasing-productivity-and-confidence-for-n-linked-glycan-analysis-of-biosimilars-using-the-bioaccord-system.html>>
6. Shion, H.; *et al.* Analysis of Antibody Drug Conjugates (ADCs) by Native Spectrometry on the BioAccord System. Waters Technology Brief. 2019. 720006570EN. <<https://www.waters.com/nextgen/us/en/library/application-notes/2019/analysis-of-antibody-drug-conjugates-adcs-by-native-mass-spectrometry-on-the-bioaccord-system.html>>
7. Ranbaduge, N.; *et al.* Establishing an Integrated Peptide Attribute Profiling and Monitoring Workflow for Improved Productivity and Confidence. Waters Application Note. 2019. 720006602EN <

<https://www.waters.com/nextgen/us/en/library/application-notes/2019/integrated-peptide-attribute-profiling-monitoring-workflow-for-improved-productivity-confidence.html> .

8. Ippoliti, S.; *et al.* Online IEX-MS of mAb Charge Variants using a BioResolve SCX mAb Column, IonHance CS-MS pH Concentrates, and BioAccord System. Waters Application Note. 2019. 720006672EN <<https://www.waters.com/nextgen/us/en/library/application-notes/2019/online-iex-ms-mab-charge-variants-bioresolve-scx-mab-column-ionhance-cx-ms-ph-concentrates-and-bioaccord-system.html>> .
9. Shion, H.; *et al.* Evaluating BioAccord LC-MS System Performance for Deployment of Compact SmartMS Enabled Biopharma Solution in Late Stage Development and QC Environment – Intact Mass Analysis. Waters Application Note. 2020. 720006722EN <<https://www.waters.com/nextgen/us/en/library/application-notes/2020/bioaccord-lc-ms-system-performance-smartms-enabled-biopharma-solution-development-qc-intact-mass-analysis.html>> .
10. Shion, H.; *et al.* Evaluating BioAccord System performance for deployment of compact smart MS system in late stage development and QC environment – peptide mapping analysis. Waters Application Note. 2020. 720006831EN <<https://www.waters.com/nextgen/us/en/library/application-notes/2020/evaluating-bioaccord-system-performance-in-late-stage-development-and-qc-environment-peptide-mapping-analysis.html>> .
11. Shion H, Birdsall R, Cubbedge S, and Chen W. Development of Integrated Informatics Workflows for the Automated Assessment of Comparability for Antibody Drug Conjugates (ADCs) Using LC-UV and LC-UV/MS. Waters Application Note. 2015 May. 720005306EN <<https://www.waters.com/nextgen/us/en/library/application-notes/2015/routine-ms-detection-for-usp-chromatographic-methods.html>> .

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© 2023 Waters Corporation. All Rights Reserved.

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