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
Report Content

Select the sections to include in your report:

- Injection List
- Blanks
- System Suitability
- Monitored Attributes
- New Peak Detection
- Acquisition Record

Generate Report

Peptide MAM Report

Generated By: waters_connect Administrator
 Generated Date: 29 Sep 2020 12:46
 Processing Date: 25 Sep 2020 13:06

Analysis
 Name: Peptide MAM
 Version: 2
 Location: Company\MAM\Demo data
 Description: bundle 34

Processing Method
 Name: Peptide MAM
 Version: 2
 Location: Company\MAM\Demo data
 Description: No Description

Injection List

| # | Injection Name | Type | Acquisition Status |
|----|-------------------------------|---------|--------------------|
| 1 | Blank_1 | Blank | Complete |
| 2 | System suitability peptides_1 | SST | Complete |
| 3 | System suitability peptides_2 | SST | Complete |
| 4 | System suitability peptides_3 | SST | Complete |
| 5 | Blank_2 | Blank | Complete |
| 6 | Control | Unknown | Complete |
| 7 | Unknown sample_1 | Unknown | Complete |
| 8 | Unknown sample_2 | Unknown | Complete |
| 9 | Unknown sample_3 | Unknown | Complete |
| 10 | Blank_3 | Blank | Complete |
| 11 | System suitability peptides_4 | SST | Complete |
| 12 | System suitability peptides_5 | SST | Complete |
| 13 | System suitability peptides_6 | SST | Complete |
| 14 | Blank_4 | Blank | Complete |

Results - % Modified

| Injection Name | WSV% (Relative to Control) | WSV% (Relative to Control) | WSV% (Relative to Control) | WSV% (Relative to Control) | WSV% (Relative to Control) |
|------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|
| Control | 1.30 | 0.91 | 0.65 | 1.47 | 41.32 |
| Unknown sample_1 | 1.18 | 1.18 | 0.68 | 4.12 | 42.29 |
| Unknown sample_2 | 1.24 | 0.90 | 0.62 | 1.30 | 42.10 |
| Unknown sample_3 | 1.35 | 1.49 | 0.66 | 4.14 | 42.44 |

New Peak Detection

| Injection Name | Relative to Control | Relative to Control | Relative to Control | Relative to Control | Relative to Control |
|------------------|---------------------|---------------------|---------------------|---------------------|---------------------|
| Control | 1.64 | 0.60 | 1.30 | | |
| Unknown sample_1 | 1.52 | 0.75 | 3.84 | | |
| Unknown sample_2 | 0.69 | 0.17 | 1.76 | | |
| Unknown sample_3 | 1.81 | 0.26 | 3.99 | | |

Processing Parameters

| Parameter | Value |
|-----------------------------|-------|
| System suitability criteria | 0.05 |
| System suitability criteria | 0.05 |

ACQUITY HDA Detector

Information
 Serial Number: 100047
 Instrument Serial Number: 2020

Instrument Setup
 Model: HDA
 Configuration: HDA-2000
 Software Version: 3.1.2
 Firmware Version: 1.7.307 (Rev. 14.03.08)
 Manufacturer: Waters Corporation
 Manufacturer Part Number: 640 p/n
 Manufacturer Part Number: 2000 p/n
 Manufacturer Part Number: 1000 p/n
 Manufacturer Part Number: 1000 p/n

Sample Manager FTM
 Software Version: 3.1.2
 Configuration: HDA-2000
 Model: FTM
 Manufacturer: Waters Corporation
 Manufacturer Part Number: 640 p/n
 Manufacturer Part Number: 2000 p/n
 Manufacturer Part Number: 1000 p/n
 Manufacturer Part Number: 1000 p/n

System suitability peptides_2

Experimental Record
 HDA Detector
 Manufacturer: Waters Corporation
 Model: HDA-2000
 Serial Number: 100047
 Instrument Serial Number: 2020
 Configuration: HDA-2000
 Software Version: 3.1.2
 Firmware Version: 1.7.307 (Rev. 14.03.08)
 Manufacturer: Waters Corporation
 Manufacturer Part Number: 640 p/n
 Manufacturer Part Number: 2000 p/n
 Manufacturer Part Number: 1000 p/n
 Manufacturer Part Number: 1000 p/n

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1. Xu, Weichen *et al.* A Quadrupole Dalton-based Multi-Attribute Method for Product Characterization, Process Development, and Quality Control of Therapeutic Proteins. *mAbs* vol. 9,7 (2017): 1186–1196. doi:10.1080/19420862.2017.1364326.
2. Alt N, Zhang TY, Motchnik P, Taticek R, Quarmby V, Schlothauer T, Beck H, Emrich T, Harris RJ. Determination of Critical Quality Attributes for Monoclonal Antibodies using Quality by Design Principles. *Biologicals*. 2016;44:291–305. doi:10.1016/j.biologicals.2016.06.005. PMID:27461239.
3. Finkler C, Krummen L. Introduction to the Application of QbD Principles for the Development of Monoclonal Antibodies. *Biologicals*. 2016;44:282-90. doi:10.1016/j.biologicals.2016.07.004. PMID:27546184.
4. Rogers, Richard S *et al.* Development of a Quantitative Mass Spectrometry Multi-Attribute Method for Characterization, Quality Control Testing and Disposition of Biologics. *mAbs* vol. 7,5 (2015): 881-90. doi:10.1080/19420862.2015.1069454.
5. Li, W, Kerwin, J. L, Schiel, J, Formolo, T, Davis, D, Mahan, A, Benchaar, S. A. State-of-the-Art and Emerging Technologies for Therapeutic Monoclonal Antibody Characterization Volume 2. *Biopharmaceutical Characterization: The NISTmAb Case Study*. January 1, 2015 , 119–183.

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[ACQUITY UPLC I-Class PLUS系统 <https://www.waters.com/134613317>](https://www.waters.com/134613317)

[ACQUITY UPLC可变波长紫外检测器 <https://www.waters.com/514228>](https://www.waters.com/514228)

[waters_connect HUB <https://www.waters.com/waters/nav.htm?cid=135040165>](https://www.waters.com/waters/nav.htm?cid=135040165)

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© 2021 Waters Corporation. All Rights Reserved.