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The screenshot displays the UNIFI software interface for the PFTreDA database. On the left, a search bar and a list of compounds are visible, with 'PFTreDA' selected. The main window shows the details for 'PFTreDA (119820 PFAS Library)'. The 'Property Value' table includes:

Property	Value
Item type	Compound
Item description	PFAS: carboxylate: legacy
IUPAC name	
Formula	C14HF27O2
Hill formula	C14HF27O2
Average molar mass	714.1134
Monoisotopic mass	713.9545
Item tag	
InChI	15/C14HF27O2/c15-2(16.1, 42)43(17.18)4(19.20)5 (21.22)6(23.24)7(25.26)8 (27.28)9(29.30)10(31.32)11 (33.34)12(35.36)13(37.38) 14(39.40)41/h(42,43)

The chemical structure of Perfluorotetradecanoic acid is shown as a ball-and-stick model. Below the structure is a 'Detection results' table:

Priority	Neutral Mass (Da)	Adduct	Charge	Fragmentation type	Expected m/z	Observed...	Expected RT (min)	Ionization technique	Detail type
Instrument model: Unknown, Instrument serial no.: Manually created, Created by administrator on Aug 11, 2020 (1 items)									
Imported from spreadsheet									
1	713.9545	-H	-1	None	712.9473		15.320	ESI-	MSe
2				CID	668.9574		15.320	ESI-	MSe
4				CID	218.9862		15.320	ESI-	MSe
3				CID	168.9894		15.320	ESI-	MSe

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Names Properties Searches Spectra Vendors Articles More ▼

Names and Synonyms

Validated by Experts, Validated by Users, Non-Validated, Removed by Users EDIT

**2,2,3,3,4,4-Hexafluor-4-(trifluoromethoxy)butansäure** [German] [ACD/IUPAC Name]  
**2,2,3,3,4,4-Hexafluoro-4-(trifluoromethoxy)butanoic acid** [ACD/IUPAC Name]  
**863090-89-5** [RN]  
**Acide 2,2,3,3,4,4-hexafluoro-4-(trifluorométhoxy)butanoïque** [French] [ACD/IUPAC Name]  
**Butanoic acid, 2,2,3,3,4,4-hexafluoro-4-(trifluoromethoxy)-** [ACD/Index Name]  
**Perfluoro(4-methoxybutanoic) acid**  
 MFCD07784238 [MDL number]  
 Perfluoro-4-methoxybutanoic acid  
 Perfluoro-4-methoxybutanoic acid (PFMOBA)

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---

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1. V Mulabagal, L Liu, J Qi, C Wilson, JS Hayworth. A Rapid UHPLC-MS/MS Method for Simultaneous Quantitation of 23 Perfluoroalkyl Substances (PFAS) in Estuarine Water. *Talanta*.2018, 190, 95–102.
2. RC Buck, J Franklin, U Berger, JM Conder, IT Cousins, P de Voogt, AA Jensen, K Kannan, SA Mabury, SPJ van Leeuwen. Perfluoroalkyl and Polyfluoroalkyl Substances in the Environment: Terminology, Classification, and Origins. *Integr Environ Assess Manag*, 2011, 7,4, 513–541.
3. Md. Al Amin, Z Sobhani, Y Liu, R Dharmaraja, S Chadalavada, Ravi Naidua, JM. Chalker, C Fang. Recent Advances in the Analysis of Per- and Polyfluoroalkyl Substances (PFAS) – a Review. *Environmental Technology and Innovation*.2020, 19, 100879.
4. S Ullah, T Alsberg, U Berger. Simultaneous Determination of Perfluoroalkyl Phosphonates, Carboxylates, and Sulfonates in Drinking Water. *J Chrom A*.2011, 1218, 6388-6395.
5. M Strynar, S Dagino, R McMahan, S Liang, A Lindstrom, E Andersen, L McMillian, M Thurman, I Ferrer, C Ball. Identification of Novel Perfluoroalkyl Carboxylic Acid (PFECAs) and Sulfonic Acids (PFESAs) in Natural Waters Using Accurate Mass Time-of-flight Mass Spectrometry (TOF MS). *Environ Sci Technol*.2015, 49, 11622-11630.
6. J McCord, M Strynar. Identification of Per- and Polyfluoroalkyl Substances in the Cape Fear River by High Resolution Mass Spectrometry and Nontargeted Screening. *Environ Sci Technol*.2019, 53, 4717-4727.
7. Y Liu, LA D' Agostino, G Qu, G Jiang, JW Martin. High Resolution Mass Spectrometry Methods for Nontarget Discovery and Characterization of Poly- and Per-fluoroalkyl Substances (PFASs) in Environmental and Human Samples. *Trends Analyt Chem*.2019, 121, 115420.

8. WA Gebbink, L van Asseldonk, SPJ van Leeuwen. Presence of Emerging Per- and Polyfluoroalkyl Substances (PFASs) in River and Drinking Water near a Fluorochemical Production Plant in the Netherlands. *Environ Sci Technol.* 2017, 51, 11057–11065.
9. S Newton, R McMahan, JA Stoeckel, M Chislock, A Lindstrom, M Strynar. Novel Polyfluorinated Compounds Identified Using High Resolution mMass Spectrometry Downstream of Manufacturing Facilities near Decatur Alabama. *Environ Sci Technol.* 2016, 51, 3, 1544–1552.
10. L Mullin, DR. Katz, N Riddell, R Plumb, JA Burgess, LWY Yeung, I Ericson Jogsten. Analysis of Hexafluoropropylene Oxide-Dimer Acid (HFPO-DA) by Liquid Chromatography-Mass Spectrometry (LC-MS): Review of Current Approaches and Environmental Levels. *Trends Analyt Chem*, 2019, 118, 828–839.
11. The Organisation for Economic Co-operation and Development (OECD). Toward a New Comprehensive Global Database of Per- and Polyfluoroalkyl Substances (PFASs): Summary Report on Updating the OECD 2007 list of Per- and Polyfluoroalkyl Substances (PFASs), OECD Ser. Risk Manag. 2018.
12. L Mullin, J Burgess. 采用PFAS分析试剂盒对全氟烷基化合物(PFAS)进行超低浓度检测. 沃特世公司技术简报. 720005701ZH <<https://www.waters.com/webassets/cms/library/docs/720005701en.pdf>> .
13. ASTM D7968, Standard Test Method for Determination of Polyfluorinated Compounds in Soil by Liquid Chromatography Tandem Mass Spectrometry (LC-MS/MS), ASTM International, West Conshohocken, PA, 2017, <https://www.astm.org> <<https://www.astm.org>> .
14. K Organtini, K Rosnack. A Large Volume Injection Technique Using Simplified Sample Preparation for Perfluorinated Alkyl Substance in Soils in Accordance with ASTM 7968. Waters Application Note 720006764EN <<https://www.waters.com/webassets/cms/library/docs/720006764en.pdf>> .
15. K Organtini, G Cleland, K Rosnack, 参照ASTM 7979-17利用大体积直接进样法分析环境水样品中的全氟烷基化合物(PFAS). 720006329ZH <<https://www.waters.com/webassets/cms/library/docs/720006329en.pdf>> .
16. An Overview of the Principles of MS<sup>E</sup>, the Engine that Drives MS Performance. Waters White Paper 720004036EN <<https://www.waters.com/webassets/cms/library/docs/720004036en.pdf>> .
17. U.S. Environmental Protection Agency. CompTox Chemicals Dashboard. [https://comptox.epa.gov/dashboard/chemical\\_lists/pfasmaster](https://comptox.epa.gov/dashboard/chemical_lists/pfasmaster) <[https://comptox.epa.gov/dashboard/chemical\\_lists/pfasmaster](https://comptox.epa.gov/dashboard/chemical_lists/pfasmaster)> . Accessed September 2020.

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© 2021 Waters Corporation. All Rights Reserved.