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
1. Higel F, et al. N-glycans of Complex Glycosylated Biopharmaceuticals and their Impact on Protein Clearance. *European Journal of Pharmaceutics and Biopharmaceutics*. 2019; 139, 123–131.
2. Zhang P, et al. Challenges of Glycosylation Analysis and Control: An Integrated Approach to Producing Optimal and Consistent Therapeutic Drugs. *Drug Discovery Today*. 2016; 21(5), 740–765.
3. Yu, Y. 将新型游离寡糖标记试剂RapiFluor-MS与集成式UPLC-FLR/QToF MS系统相结合用于分析低丰度N-糖. 沃特世应用纪要, 2015. [720005383ZH < https://www.waters.com/nextgen/us/en/library/application-notes/2015/applying-glycan-tagging-reagent-rapi-fluor-ms-and-integrated-uplc-flr-qtof-ms-system-for-low-abundant-n-glycan-analysis.html > .](https://www.waters.com/nextgen/us/en/library/application-notes/2015/applying-glycan-tagging-reagent-rapi-fluor-ms-and-integrated-uplc-flr-qtof-ms-system-for-low-abundant-n-glycan-analysis.html)
4. Zhang X, Reed C, Shion H, Alley W, Birdsall R, Yu Y. 使用BioAccord系统提高生物类似药N-糖分析的工作效率和可信度. 沃特世应用纪要, 2019. [720006545ZH < https://www.waters.com/nextgen/us/en/library/application-notes/2019/increasing-productivity-and-confidence-for-n-linked-glycan-analysis-of-biosimilars-using-the-bioaccord-system.html > .](https://www.waters.com/nextgen/us/en/library/application-notes/2019/increasing-productivity-and-confidence-for-n-linked-glycan-analysis-of-biosimilars-using-the-bioaccord-system.html)
5. DeLoff M. A Comparison of Premier Columns with MaxPeak HPS Technology Versus PEEK-lined Column Hardware. *Waters Application Note*, 2021. [720007210EN < https://www.waters.com/nextgen/us/en/library/application-notes/2021/a-comparison-of-premier-columns-with-maxpeak-hps-technology-versus-peek-lined-column-hardware.html > .](https://www.waters.com/nextgen/us/en/library/application-notes/2021/a-comparison-of-premier-columns-with-maxpeak-hps-technology-versus-peek-lined-column-hardware.html)

6. Isaac G., Plumb R. ACQUITY Premier液相色谱技术显著改善磷酸化和羧化脂质的灵敏度、峰形和回收率.沃特世应用纪要, 2021.720007092ZH <<https://www.waters.com/nextgen/us/en/library/application-notes/2021/acquity-premier-lc-technology-significantly-improves-sensitivity,-peak-shape-and-recovery-for-phosphorylated-and-carboxylate-lipids.html>>
7. Tanna N. and Plumb R. 采用配备MaxPeak HPS技术的ACQUITY Premier系统和色谱柱进行吉非替尼（一种EGFR抑制剂）生物分析的优势.沃特世应用纪要, 2021.720007122ZH <<https://www.waters.com/nextgen/us/en/library/application-notes/2021/advantages-of-using-acquity-premier-system-and-columns-with-maxpeak-hps-technology-for-bioanalysis-of-gefitinib-an-egfr-inhibitor.html>>
8. Liu X. and Lauber M. High Performance Surface Technology Improves HILIC Profiling of Released N-Glycans.Waters Application Note, 2021.
9. Kang J.F., et al. Lysosomal Targeting Enhancement by Conjugation of Glycopeptides Containing Mannose-6-phosphate Glycans Derived from Glyco-engineered Yeast.Scientific Reports. 2018; 8, 8730.
10. Yu Y. N-糖的表征与分析：将精确质量、标准GU值和功能强大的UNIFI软件有机结合，更准确地进行寡糖结构确证.沃特世应用纪要, 2016.720004845ZH <<https://www.waters.com/nextgen/us/en/library/application-notes/2016/n-linked-glycan-accurate-mass-reference-glucose-units-glycan-assignments0.html>> .

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© 2021 Waters Corporation. All Rights Reserved.