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---

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1. European Monitoring Centre for Drugs and Drug Addiction (2020). *New psychoactive substances: global markets, global threats and COVID-19 pandemic. An update from the EU Early warning System* (December 2020), Publications of the European Union, Luxembourg.  
[https://www.emcdda.europa.eu/publications/rapid-communication/new-psychoactive-substances-global-markets-global-threats-and-covid-19-pandemic\\_en](https://www.emcdda.europa.eu/publications/rapid-communication/new-psychoactive-substances-global-markets-global-threats-and-covid-19-pandemic_en) <  
[https://www.emcdda.europa.eu/publications/rapid-communication/new-psychoactive-substances-global-markets-global-threats-and-covid-19-pandemic\\_en](https://www.emcdda.europa.eu/publications/rapid-communication/new-psychoactive-substances-global-markets-global-threats-and-covid-19-pandemic_en)> (accessed 24 Aug 2022).
2. Wood M. 搭载LiveID的RADIANT ASAP – 实现快速、特异、简便的药物筛查. 沃特世应用纪要. 720007125ZH, 2021.
3. Evans-Brown M., *et al.* EMCDDA Technical Report on the New Psychoactive Substance *N,N*-diethyl-2-[[4-(1-methylethoxy)phenyl]methyl]-5-nitro-1*H*-benzimidazole-1-ethanamine (isotonitazene).  
[https://www.emcdda.europa.eu/publications/initial-reports/isotonitazene\\_en](https://www.emcdda.europa.eu/publications/initial-reports/isotonitazene_en) <  
[https://www.emcdda.europa.eu/publications/initial-reports/isotonitazene\\_en](https://www.emcdda.europa.eu/publications/initial-reports/isotonitazene_en)> (accessed 24 Aug 2022).
4. Advisory Council on the Misuse of Drugs. ACMD report – A Review of the Evidence on the Use and harm of 2-Benzyl Benzimidazole (‘nitazene’) and Piperidine Benzimidazolone (‘bromophine-like’) Opioids. July 2022. <https://www.gov.uk/government/publications/acmd-advice-on-2-benzyl-benzimidazole-and-piperidine-benzimidazolone-opioids/acmd-advice-on-2-benzyl-benzimidazole-and-piperidine-benzimidazolone-opioids-accessible-version> <  
<https://www.gov.uk/government/publications/acmd-advice-on-2-benzyl-benzimidazole-and-piperidine-benzimidazolone-opioids/acmd-advice-on-2-benzyl-benzimidazole-and-piperidine-benzimidazolone-opioids-accessible-version>>

[piperidine-benzimidazolone-oids/acmd-advice-on-2-benzyl-benzimidazole-and-piperidine-benzimidazolone-oids-accessible-version](#)> (accessed 24 Aug 2022).

5. Drug Enforcement Administration, Diversion Control Division, Drug and Chemical Evaluation Section, May 2022. [https://www.deadiversion.usdoj.gov/drug\\_chem\\_info/benzimidazole-oids.pdf](https://www.deadiversion.usdoj.gov/drug_chem_info/benzimidazole-oids.pdf) <[https://www.deadiversion.usdoj.gov/drug\\_chem\\_info/benzimidazole-oids.pdf](https://www.deadiversion.usdoj.gov/drug_chem_info/benzimidazole-oids.pdf)> (accessed 24 Aug 2022).
6. Canadian Centre on Substance Use and Addiction. CCENDU Drug Alert Nitazenes. March 2022. [https://www.ccsa.ca/sites/default/files/2022-03/CCSA-CCENDU-Drug-Alert-Nitazenes-2022-en\\_0.pdf](https://www.ccsa.ca/sites/default/files/2022-03/CCSA-CCENDU-Drug-Alert-Nitazenes-2022-en_0.pdf) <[https://www.ccsa.ca/sites/default/files/2022-03/CCSA-CCENDU-Drug-Alert-Nitazenes-2022-en\\_0.pdf](https://www.ccsa.ca/sites/default/files/2022-03/CCSA-CCENDU-Drug-Alert-Nitazenes-2022-en_0.pdf)> (accessed 24 Aug 2022).
7. Di Trana A, et al. Synthetic Benzimidazole Opioids: The Emerging Health Challenge for European Drug Users. *Frontiers in Psychiatry*. 13:858234. 2022.
8. Lee E, Cooper J, Martin N, Wood M. 更新RADIANT ASAP查获药物参比数据库, 沃特世白皮书, 720007838ZH <<https://www.waters.com/webassets/cms/library/docs/720007838en.pdf>> 2023.
9. Mistry N.S, Cooper J, Wood M. MRM毒理学筛查方法配合Waters Xevo TQ-S micro的运用扩展. 沃特世应用纪要. 720007748ZH. 2022年.

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© 2023 Waters Corporation. All Rights Reserved.