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---

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1. Schymanski EL, Jeon J, Gulde R, Fenner K, Ruff M, Singer HP, Hollender J. Identifying Small Molecules via High Resolution Mass Spectrometry: Communicating Confidence. *Environ.Sci.Technol.* 2014, 48 (4), 2097–2098.
2. Cooper-Shepherd DA, Wildgoose J, Kozlov B, Johnson WJ, Tyldesley-Worster R, Palmer ME, Hoyes JB, McCullagh M, Jones E, Tonge R, Marsden-Edwards E, Nixon P, Verenchikov A, Langridge JI. Novel Hybrid Quadrupole-Multireflecting Time-of-Flight Mass Spectrometry System. *J Am Soc Mass Spectrom* .2023 Feb 1;34(2):264–272.
3. McCullagh M, Palmer M, Marsden-Edwards E, Cooper-Shepherd D和James Langridge, 使用配备UPLC™的多反射飞行时间质谱仪探索ppb质量精度对代谢物鉴定的影响（A部分），沃特世应用纪要，[720007896ZH](#)，2023年4月。
4. FDA, US. "Guidance for industry: safety testing of drug metabolites." *Center for Drug Evaluation and Research, Rockville, MD, USA* (2008).

5. Bradshaw PR, Athersuch TJ, Stachulski AV, Wilson ID. Acyl glucuronide reactivity in perspective, *Drug Discovery Today*, Volume 25, Issue 9, 2020, Pages 1639–1650.
6. Mortensen RW, Corcoran O, Cornett C, Sidelmann UG, Troke J, Lindon JC, Nicholson JK, Hansen SH. LC-1H NMR used for determination of the elution order of S-naproxen glucuronide isomers in two isocratic reversed-phase LC-systems. *J Pharm Biomed Anal.* 2001 Jan;24(3):477–85.
7. Higton D, Wilson ID, Vissers JPC, Plumb RS. An automated method to study the rapid intramolecular transacylation of drug acyl glucuronides using Cyclic Ion Mobility Spectrometry-Mass Spectrometry Waters Corporation Poster: [PSTR135089718 <https://www.waters.com/webassets/cms/library/docs/2021asms\\_higton\\_acylglucuronide.pdf>](https://www.waters.com/webassets/cms/library/docs/2021asms_higton_acylglucuronide.pdf).
8. Falany CN, Ström P, Swedmark S. Sulphation of o-desmethylnaproxen and related compounds by human cytosolic sulfotransferases. *Br J Clin Pharmacol.* 2005 Dec;60(6):632–40.
9. Kiang CH, Lee C, Kushinsky S. Isolation and identification of 6-desmethylnaproxen sulfate as a new metabolite of naproxen in human plasma. *Drug Metab Dispos* 1989; 17: 43–8.
10. Nagao T, Yukihiro D, Fujimura Y, et al. Power of isotopic fine structure for unambiguous determination of metabolite elemental compositions: in silico evaluation and metabolomic application. *Anal Chim Acta.* 2014;813:70–76.
11. Ipsen A. Efficient calculation of exact fine structure isotope patterns via the multidimensional Fourier transform. *Anal Chem.* 2014;86(11):5316–5322.
12. Lertratanakoon K and Horning M G. Metabolism of carbamazepine. *Drug Metabolism and Disposition* January 1982, 10 (1) 1–10.
13. Bahlmann A, Brack W, Schneider RJ, Krauss M. Carbamazepine and its metabolites in wastewater: Analytical pitfalls and occurrence in Germany and Portugal. *Water Res.* 2014 Jun 15;57:104–14. doi: 10.1016/j.watres.2014.03.022. Epub 2014 Mar 19. PMID: 24704908.
14. Richter WJ, Kriemler P, Faigle JW. (1978). Newer Aspects of the Biotransformation of Carbamazepine: Structural Characterization of Highly Polar Metabolites. In: Frigerio, A. (eds) *Recent Developments in Mass Spectrometry in Biochemistry and Medicine*. Springer, Boston, MA. [https://doi.org/10.1007/978-1-4613-3991-5\\_1](https://doi.org/10.1007/978-1-4613-3991-5_1)

15. Russell JL, Spiller HA, Baker DD. Markedly elevated carbamazepine-10,11-epoxide/carbamazepine ratio in a fatal carbamazepine ingestion, *Case Rep. Med.* 2015 (2015), 369707.
16. Kerr BM, Thummel KE, Wurden CJ, Klein SM, Kroetz DL, Gonzalez FJ, Levy R, Human liver carbamazepine metabolism: Role of CYP3A4 and CYP2C8 in 10,11-epoxide formation, *Biochemical Pharmacology*, Volume 47, Issue 11, 1994, Pages 1969–1979, ISSN 0006–2952.
17. Batucan NSP, Tremblay LA, Northcott GL, Matthaei CD. Medicating the environment? A critical review on the risks of carbamazepine, diclofenac and ibuprofen to aquatic organisms, *Environmental Advances*, Volume 7, 2022, 100164.
18. Gago-Ferrero p, Bletsou AA, Damalas DE, Aalizadeh R, Alygizakis NA, Singer HP, Hollender J, Thomaidis NS. Wide-scope target screening of >2000 emerging contaminants in wastewater samples with UPLC-Q-ToF-HRMS/MS and smart evaluation of its performance through the validation of 195 selected representative analytes. *Journal of Hazardous Materials*, Volume 387, 2020, 121712.
19. Carere M, Polesello S, Kase R, Gawlik BM. (2015). The Emerging Contaminants in the Context of the EU Water Framework Directive. In: Petrovic M, Sabater S, Elosegi A, Barceló, D. (eds) *Emerging Contaminants in River Ecosystems. The Handbook of Environmental Chemistry*, vol 46. *Springer, Cham*.

---

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[ACQUITY UPLC I-Class PLUS系统 <https://www.waters.com/134613317>](https://www.waters.com/134613317)

[SELECT SERIES MRT <https://www.waters.com/waters/nav.htm?cid=135082877>](https://www.waters.com/waters/nav.htm?cid=135082877)

[MassLynx MS软件 <https://www.waters.com/513662>](https://www.waters.com/513662)

[waters\\_connect <https://www.waters.com/waters/nav.htm?cid=135040165>](https://www.waters.com/waters/nav.htm?cid=135040165)

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© 2023 Waters Corporation. All Rights Reserved.

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