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The screenshot displays the INTACT Mass software interface. On the left, a sidebar contains navigation options: 'Define LC peak detection parameters', 'Define peak deconvolution parameters', 'Select modifiers' (highlighted), 'Specify acceptable limits for quantitation', and 'Specify identification and quantitation limits'. The main area is titled 'Select modifiers' and includes a search bar and a 'Save' button. Below this is a table of 'Available modifiers' with columns for Name, Average mass, Monoisotopic mass, Formula, and Type. The 'n-A' modifier is selected. To the right, the 'Selected modifiers' section shows a list of modifiers: 'n-A' and 'n+A', both with a maximum of 10. Below that, the 'Preview of possible modifications: 20' section shows a table of 20 possible modifications, including 'n-A', 'n+A', 'n-A(2)', 'n-A(3)', and 'n-A(4)'. A system message 'Activate Windows Go to Settings to activate Windows.' is visible in the bottom right of the preview table.

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1. Jackson NAC, Kester KE, Casimiro D, Gurunathan S, DeRosa F. The Promise of mRNA Vaccines: A Biotech and Industrial Perspective, *npj Vaccines*, 2020, 11, 1–6.
2. Sahin U, Kariko K, Tureci O. mRNA-Based Therapeutics – Developing a New Class of Drugs, *Nat Rev*, 2014, 13, 759–780.
3. Sheets MD, Wickens W. Two Phases in the Addition of a Poly(A) Tail, *Genes & Develop*, 1989, 3, 1401–1412.
4. Dreyfus M, Regnier P. The Poly(A) Tail of mRNAs: Bodyguard in Eukaryotes, Scavenger in Bacteria, *Cell*, 2002, 111, 611–613.
5. Chang H, Lim J, Ha M, Kim VN. TAIL-seq: Genome-Wide Determination of Poly(A) Tail Length and 3' End Modifications, *Mol Cell*, 2014, 53, 1044–1052.
6. Nicholson AL, Pasquinelly AE. Tales of Detailed Poly(A) Tails, *Trends Cell Biol*, 2019, 29, 191–200.
7. Jiang T, Yu N, Kim J, Murgu JR, Kissai M, Ravichandran K, Miracco EJ, Presnyak V, Hua S. Oligonucleotide Sequence Mapping of Large Therapeutic mRNAs via Parallel Ribonuclease Digestions and LC-MS/MS, *Anal Chem*, 2019, 91, 8500–8506.
8. Vanhinsbergh C, Criscuolo A, Sutton JN, Murphy K, Williamson AJK, Cook K, Dickman MJ. Characterization and Sequence Mapping of Large RNA and mRNA Therapeutics Using Mass Spectrometry, *Anal Chem*, 2022, 94, 7339–7349.
9. Muttach F, Muthmann N, Rentmeister A. Synthetic mRNA Capping, *Beilstein J Org Chem*, 2017, 13, 2819–2832.
10. Beverly M, Dell A, Parmar P, Houghton L. Label-Free Analysis of mRNA Capping Efficiency Using RNase H probes and LC-MS, *Anal Bioanal Chem*, 2016, 408, 5021–5030.
11. Nguyen J, Chan SH, Roy B, Gilar M, Robb B, Chen W, Lauber MA. Rapid Analysis of Synthetic mRNA Cap Structure Using Ion-Pairing RPLC with the BioAccord LC-MS System, Waters Application Note, [720007329](#), 2021.
12. Muthmann N, Spacek P, Reichert D, van Dulmen M, Rentmeister A. Quantification of mRNA Cap-Modifications by Means of LC-QqQ-MS, *Methods*, 2022, 203, 196–206.
13. Chan SH, Whipple JM, Dai N, Kelley TM, Withers K, Tzertzinis G, Correa Jr IR, Robb GB. RNase-H-based Analysis of Synthetic mRNA 5' Cap Incorporation, *RNA*, 2022, 28, 1144–1155.

14. Beverly M, Hagen C, Slack O. Poly A Tail Length Analysis of In Vitro Transcribed mRNA by LC-MS, *Anal Bioanal Chem*, 2018, 410, 1667–1677.
15. Strezsak SR, Pimentel AJ, Hill IT, Beuning PJ, Skizim NJ. Novel Mobile Phase to Control Charge States and Metal Adducts in the LC/MS for mRNA Characterization Assays, *ACS Omega Reviews*, 2022, 7, 22181–22191.
16. Doneanu C E, Fox J, Harry E, Knowles C, Yu Y Q, Fredette J, Chen W. An Automated Compliance-Ready LC-MS Workflow for Intact Mass Confirmation and Purity Analysis of Oligonucleotides, Waters Application Note, [720006820](#), 2020.
17. Doneanu C E, Fox J, Harry E, Knowles C, Yu Y Q, Fredette J, Chen W. Intact Mass Confirmation Analysis on the BioAccord LC-MS System for a Variety of Extensively Modified Oligonucleotides, Waters Application Note, [720007028](#), 2020.
18. Doneanu C E, Knowles C, Fox J, Harry E, Yu Y Q, Fredette J, Chen W. Analysis of Oligonucleotide Impurities on the BioAccord System with ACQUITY Premier, Waters Application Note, [720007301](#), 2021.
19. Doneanu C E, Boyce P, Shion H, Fredette J, Berger S J, Gastall H, Yu Y Q. LC-MS Analysis of siRNA, Single Guide RNA and Impurities Using the BioAccord™ System with ACQUITY™ Premier and New Automated INTACT Mass Application, Waters Application Note, [720007546](#), 2022.
20. Gilar M. Size-Exclusion Chromatography Method for Poly(A) Tail Analysis of mRNA, Waters Application Note, [720007853](#), 2023.
21. Shion H, Boyce P, Berger S J, Yu Y Q. INTACT Mass™ - a Versatile waters_connect™ Application for Rapid Mass Confirmation and Purity Assessment of Biotherapeutics, Waters Application Note, [720007547](#), 2022.
22. Ferrige AG, Seddon MJ, Green BN, Jarvis SA, Skilling J. Disentangling Electrospray Spectra With Maximum Entropy, *Rapid Commun. Mass Spectrom.*, 1992, 6, 707–711.

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© 2023 Waters Corporation. All Rights Reserved.

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