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1. Cooper-Shepherd D, Wildgoose J, Kozlov B, Johnson WJ, Tyldesley-Worster R, Palmer M, Hoyes J, McCullagh M, Jones E, Tonge R, Marsden-Edwards E, Nixon P, Verenchikov A, and Langridge J. Novel Hybrid Quadrupole-Multireflecting Time-of-Flight Mass Spectrometry System. *Journal of the American Society for Mass Spectrometry* 2023 34 (2), 264–272.
2. McCullagh M, Palmer M, Eatough D, Marsden-Edwards E, Langridge J, Cooper-Shepherd D. 使用配备 UPLC™的多反射飞行时间质谱仪探索ppb质量精度对代谢物鉴定的影响 (A部分) . [720007896ZH](#), 2023年4月.
3. McCullagh M, Palmer M, Eatough D, Marsden-Edwards E, Langridge J, Cooper-Shepherd D. Exploring the Impact of Part Per Billion Mass Accuracy for Metabolite Identification using Multi Reflecting Time-Of-Flight MS with UPLC Part B. [720007889](#), April 2023.
4. Falany CN, Ström P, Swedmark S. Sulphation of o-desmethylnaproxen and related compounds by human cytosolic sulfotransferases. *Br J Clin Pharmacol.* 2005 Dec;60(6):632–40.
5. Kiang CH, Lee C, Kushinsky S. Isolation and identification of 6- desmethylnaproxen sulfate as a new metabolite of naproxen in human plasma. *Drug Metab Dispos* 1989; 17: 43–8.
6. Nagao T, Yukihiro D, Fujimura Y, et al. Power of isotopic fine structure for unambiguous determination of metabolite elemental compositions: in silico evaluation and metabolomic application. *Anal Chim Acta.* 2014;813:70–76.
7. Lertratanakoon K and Horning M G. Metabolism of carbamazepine. *Drug Metabolism and Disposition* January 1982, 10 (1) 1–10.
8. Bahlmann A, Brack W, Schneider RJ, Krauss M. Carbamazepine and its metabolites in wastewater: Analytical pitfalls and occurrence in Germany and Portugal. *Water Res.* 2014 Jun 15;57:104–14. doi: 10.1016/j.watres.2014.03.022. Epub 2014 Mar 19. PMID: 24704908.
9. Richter WJ, Kriemler P, Faigle JW. (1978). Newer Aspects of the Biotransformation of Carbamazepine: Structural Characterization of Highly Polar Metabolites. In: *Frigerio, A. (eds) Recent Developments in Mass Spectrometry in Biochemistry and Medicine.* Springer, Boston, MA. https://doi.org/10.1007/978-1-4613-3991-5_1
10. Russell JL, Spiller HA, Baker DD. Markedly elevated carbamazepine-10,11-epoxide/carbamazepine ratio in a fatal carbamazepine ingestion, Case Rep. *Med.* 2015 (2015), 369707.

11. Kerr BM, Thummel KE, Wurden CJ, Klein SM, Kroetz DL, Gonzalez FJ, Levy R, Human liver carbamazepine metabolism: Role of CYP3A4 and CYP2C8 in 10,11-epoxide formation, *Biochemical Pharmacology*, Volume 47, Issue 11, 1994, Pages 1969-1979, ISSN 0006-2952.
12. Loos, M., Gerber, C., Corona, F., Hollender, J., Singer, H. (2015). Accelerated isotope fine structure calculation using pruned transition trees, *Analytical Chemistry* 87(11), 5738-5744.

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[ACQUITY UPLC I-Class PLUS系统 <https://www.waters.com/134613317>](https://www.waters.com/134613317)

[SELECT SERIES MRT <https://www.waters.com/waters/nav.htm?cid=135082877>](https://www.waters.com/waters/nav.htm?cid=135082877)

[MassLynx MS软件 <https://www.waters.com/513662>](https://www.waters.com/513662)

[waters_connect <https://www.waters.com/waters/nav.htm?cid=135040165>](https://www.waters.com/waters/nav.htm?cid=135040165)

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© 2024 Waters Corporation. All Rights Reserved.

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